

Technical data sheet

CAD/CAM Restoratives **VITA ENAMIC®**



Manufacturer

**VITA Zahnfabrik H. Rauter
GmbH & Co.KG**
P.O. Box 1338
79704 Bad Säckingen
Germany

VITA Zahnfabrik is
certified according
to:

- DIN EN ISO 13485
- RL 93/42/ECC (CE 0124)
- ISO 9001

■ Description

CAD/CAM restoration material **VITA ENAMIC®**

Hybrid ceramic blocks

■ Description

VITA ENAMIC® is a hybrid ceramic with a dual network structure that combines the best characteristics of ceramic and composite. In this dental material, the dominant fine-structure ceramic network (~ 86% by weight) is reinforced by an acrylate polymer network (~14% by weight), with both networks penetrating completely. This innovative hybrid ceramic guarantees a unique balance between high load-bearing capacity and elasticity with a high absorption potential of the chewing forces.

■ Indications

- Single tooth restorations (Inlays, onlays, veneers and crowns)
- Crown restorations in areas with high chewing force (molar range)
- Least invasive restorations
- Cosmetic veneer restorations

■ Contraindications

- Para function (bruxism)
- Bridges

■ Product properties and benefits

- Significantly lower brittleness than pure ceramics and better abrasion behaviour than composite.
- Compared to silicate ceramics, restorations require thinner wall thicknesses.
- Enamel-like antagonist-friendly abrasion properties due to the fine-structure ceramic network.
- Significantly higher elasticity than traditional dental ceramic due to the flexibility provided by the acrylate polymer network. ■ Compared to silicate ceramics, the tool life of the grinding tools is increased approx. 4-5 times.

Technical data sheet

CAD/CAM Restoratives

VITA ENAMIC®













Manufacturer

VITA Zahnfabrik H. Rauter
GmbH & Co.KG
P.O. Box 1338
79704 Bad Säckingen
Germany

VITA Zahnfabrik is
certified according to:
■ DIN EN ISO 13485
■ RL 93/42/ECC (CE 0124)
■ ISO 9001

Options

VITA ENAMIC® is offered by the manufacturer in the translucency classes HT (high translucent) and T (translucent) and offered in each case with 5 VITA SYSTEM 3D-MASTER® shades. In addition, the following table provides information about the available VITA classical shades.

0M1 = Bleach	1M1 = B1	1M2 = A1	2M2 = A2	3M2 = A3
0M1-HT 	1M1-HT 	1M2-HT 	2M2-HT 	3M2-HT 
0M1-T 	1M1-T 	1M2-T 	2M2-T 	3M2-T 

Preparation/strength of ceramic layers

■ Inlays and onlays strength of ceramic layers:

Inlays: Fissure root = at least 1.0 mm/isthmus range = at least 1.5

mm Onlays: Fissure root = at least 1.0 mm/tubercle range of the tooth= at least 1.5 mm

■ Veneers ceramic veneer thicknesses: labial: at least 0.3 mm on average/incisal third: at least 0.3 mm thicknesses of ceramic layer: average third: at least 0.3 mm/cervical third: at least 0.2 mm

■ Crowns Strength of ceramic layers:

Anterior crowns: incisal: at least 1.5 mm/circular: at least 0.8 mm

Posterior crowns: tubercle range of the tooth: at least 1.5 mm/fissure root: at least 1.0 mm circular: 0.8-1.5 mm

Technical data sheet

CAD/CAM Restoratives VITA ENAMIC®




Manufacturer

VITA Zahnfabrik H. Rauter
GmbH & Co.KG
P.O. Box 1338
79704 Bad Säckingen
Germany

VITA Zahnfabrik
is certified according to
■ DIN EN ISO 13485
■ RL 93/42/ECC (CE 0124)
■ ISO 9001

Step-by-step manufacturing process

Key
— = Praxis
■ = laboratory
■ = CADstar



Finalising/fitting of the restoration

CADstar will provide the finished and sanded restoration, already cemented; however, on request, you may do this procedure in-house.

Adhesive Attaching the restoration

- Restorations made from VITA ENAMIC® must be bonded using light-curing or dual-curing fine-hybrid composites.
- The self-adhesive composite Re-lyX Unicem (3M ESPE) is only suitable for the attachment of crowns (denture adhesion). The restoration is then treated with VITA CERAMICS ETCH for 60 seconds and subsequently silanised.
- Preferably, crowns should be bonded with a more fluid composite with dual curing character (depending on the layer thickness).
- When using more solid composite materials, the ultrasonic insertion method or a pre-heated composite can be used.
- If using thin veneers, the use of dual-curing composites should be omitted, since these can cause a slight change in colour (yellow shade) after the curing process. Therefore, a pure light-curing composite is preferable.

Cements: Fine hybrid composite with a bonding system, e.g., VITA DUO CEMENT with VITA A.R.T. BOND or PANAVIA F2.0 with ED Primer II

Technical data sheet

CAD/CAM Restoratives

VITA ENAMIC®



Manufacturer

VITA Zahnfabrik H. Rauter
 GmbH & Co.KG
 P.O. Box 1338
 79704 Bad Säckingen
 Germany

VITA Zahnfabrik is certified according to:

- DIN EN ISO 13485
- RL 93/42/ECC (CE 0124)
- ISO 9001

Material composition of the ceramic content (86% by weight or 75% by volume)

SiO ₂ (in %)	Al ₂ O ₃ (in %)	Na ₂ O (in %)	K ₂ O (in %)	B ₂ O ₃ (in %)	ZrO ₂ (in %)	KaO (in %)
58 - 63	20 - 23	9 - 11	4 - 6	0.5 - 2	< 1	< 1

Material composition of the polymer content (14% by weight or 25% by volume)

UDMA (Urethane dimethacrylate)	TEGDMA (Triethylene glycol dimethacrylate)

Physical properties (guidelines)

Density ρ (at 20°C)	2.1	[g/cm ³]
Flex resistance β (according to Schwickerath / ISO 6872)	150 - 160	[MPa] or [N/mm ²]
Elasticity modulus	30,000	[MPa] or [N/mm ²]
Abrasion	28 - 46	[μ m]
Elongation at break A	0.5	[%]
Weibull modulus	20	
Vickers hardness (HV 30)	2,500	[MPa] or [N/mm ²]
Fracture toughness K_{IC}	1.5	[MPa \sqrt{m}]