

Technical data sheet

PMMA Disc TELIO CAD®



Manufacturer

Ivoclar Vivadent AG

Bendererstrasse 2
9494 Schaan
Liechtenstein

Ivoclar Vivadent AG
is certified according
to

■ DIN EN ISO 13485
■ RL 93/42/ECC (CE 0123)

■ Designation

TELIO CAD® PMMA Disc for CAD/CAM

■ Description

TELIO CAD® discs are blanks of PMMA (polymethyl methacrylate). CAD/CAM can be used to machine single teeth as well as single or multi-unit, fully-anatomic temporary restorations. By using additional layers and stains, aesthetic optimisation can be obtained.

The most important benefits of TELIO CAD® discs are:

- high material homogeneity due to the industrial manufacturing process ■
- no toxins and benzene peroxide free
- lasting shade stability and natural fluorescence ■
- simple reproducibility of the temporary dental prosthesis

■ Indication

- Anterior and posterior crowns can be worn up to max. 12 months
- Anterior and posterior bridges with up to 2 pontics can be worn up to max. 12 months
- Implanted temporary restorations can be worn up to max. 12 months
- Therapeutic care for correction of temporomandibular joint problems and occlusal plane

■ Contraindication


- Use for definitive restorations
- Bridge constructions with more than two pontics in one piece
- Patients with parafunctions, e.g., bruxism
- if the allergy to one of the constituents is known, TELIO CAD® must not be used

■ Options

The TELIO CAD® PMMA discs are available in the following shades A1, A2, A3, A3,5 and B1.

Technical data sheet

PMMA Disc TELIO CAD®



Manufacturer

Ivoclar Vivadent AG

Bendererstrasse 2
9494 Schaan
Liechtenstein

Ivoclar Vivadent AG is certified
according to:

- DIN EN ISO 13485
- RL 93/42/ECC (CE 0123)

■ Modelling

During the construction of crown and bridge framework, the following parameters must not be their values specified: ■ Minimum wall thickness (circular): 0.8 mm

- Minimum wall thickness (occlusal): 1.5 mm
- Connector cross-section of anterior bridges with 1 pontic: min. 12 mm²
- Connector cross section of anterior bridges with 2 pontics: min 12 mm² ■
- Connector cross section of posterior bridges with 1 pontic: min. 12 mm²
- Connector cross section of posterior bridges with 2 pontics: min. 16 mm²

■ CAD/CAM Fully anatomic production with final polishing

In this processing technique, the restoration is polished and integrated into the CAD/CAM system immediately after machining. The surface finish gloss is adjusted by manual polishing. This processing step is very efficient, and makes it easy to obtain an aesthetic result quickly.

■ Finishing

The following procedure is recommended for the preparation and revision of TELIO CAD® restorations:

- Grinding the tapping point with fine cross-toothed carbide cutters
- shape corrections with fine cross-toothed carbide cutters or commercially available diamond tools
- Avoid overheating the material
- smooth the entire occlusal surface with a fine diamond tool. This will even the surface caused by CAD/CAM processing
- Observe the minimum wall thickness specifications
- Thoroughly clean the restoration before further processing (id abrasive residues remains → veneer problems!)
- Try-in and adjustment of occlusion or articulation, if necessary
- Polish (see manufacturer's instructions)

For further processing techniques, such as "cut-back" technology, see the manufacturer's processing instructions.

■ Processing instructions


Operating and processing instructions for the dental technician *IVOCLAR VIVADENT TELIO CAD® / TELIO*

CS® http://www.ivoclarvivadent.com/de/alle/produkte/chairside-cad-_-cam-bloecke/telio-cad

On the right-hand side of the processing instructions *click TELIO Lab - CAD* → Open PDF File

Technical data sheet

PMMA Disc TELIO CAD®



Manufacturer

Ivoclar Vivadent AG
Bendererstrasse 2
9494 Schaan
Liechtenstein

Ivoclar Vivadent AG is certified according:

■ **DIN EN ISO 13485**

■ **RL 93/42/ECC (CE 0123)**

■ Processing instructions

Operating and processing instructions for the dental technician *IVOCLAR VIVADENT TELIO CAD® / TELIO CS®* http://www.ivoclarvivadent.com/de/alle/produkte/chairside-cad-_-cam-bloecke/telio-cad

On the right-hand side of the processing instructions *click TELIO Lab - CAD* → Open PDF File

■ Chemical composition

PMMA (polymethyl methacrylate) 99.5% / pigments < 1.0%

■ Physical/mechanical properties (guidelines)

Density ρ (at 20°C)	1.18	[g/cm ³]
Elasticity module (at 20°C)	3,200 (± 300)	[MPa] or [N/mm ²]
Water absorption W	< 28	[µg/mm ³]
Water absorption L	< 0.6	[µg/mm ³]
Ball indentation hardness H_K	180 (± 5)	[MPa] or [N/mm ²]
Vickers hardness HV 10	190 (± 5)	[MPa] or [N/mm ²]
Flexural strength β_B	130 (± 10)	[MPa] or [N/mm ²]

■ Thermal properties (guidelines)

Vicat softening temperature (melting point) T_v	approx. 100 [°C] or 212 [°F]
Flash point T_F	> 250 [°C] or 482 [°F]
Thermal conductivity λ (at 23°C)	0.19 [W/(K*m)]
Specific heat capacity c (at 23°C)	1.47 [J/(g*K)]