Titan-Klebebasis

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Titanium bonding base

for individual zirconium oxide abutments

11.1 - series (ZTSV)

compatible with **Zimmer Dental Tapered Screw-Vent®** implants
compatible with **MIS Seven®** and **Bio Horizon®** (internal) implants

						R_x only
Instruction manual – read carefully!	Contents non sterile!	Do not reuse!	Lot number	Order number	Manufacturer	CAUTION! – Federal law (USA) restricts this product to sale by or on the order of a dentist or a physician!

Product description: The titanium base is delivered non-sterile with an abutment screw. The titanium base can be scanned with current dental scanners, and thus obtained digitally acquired geometry can be used in the manufacture of individual abutments, using CAD/CAM technologies.

Indications: The titanium base is used for customized implant-prosthetic zirconium oxide abutment manufacturing. The titanium base, after cementing with a CAD/CAM milled assembly, is the connecting component to the implant. Follow the notes in the IFU (instructions for use) from Zimmer Dental GmbH, concerning indications and contraindications.

Order No. (REF): GA_E_11-1-Serie (www.cadstar.dental/Help Center/Gebrauchsanweisungen)

Contraindications: The titanium bases of the 11.1-series can only be combined with Zimmer Dental Tapered Screw-Vent® implants. They cannot be combined with implants of a different implant type or manufacturer. The diameter of the abutment must correspond in size to the used implant in order to prevent a tissue irritation. The ti-bases and abutment screws are indicated for single use only. For fixation of the abutment on the implant, the *correct torque force of 30 Ncm*, recommended by the implant manufacturer, has to be considered carefully. Mechanical treatment of the connection part of the titanium base will damage the correct fitting of the abutment on the implant.

Composition: Titanium base, abutment screw and scan abutment: Ti6Al4V ELI ISO 5832-2 / ASTM F136

Handling method:

Ceramic abutments: Milling with CAD/CAM-machines, e.g. of zirconium oxide ceramics, according to the anatomic form of a crown or coping. The ceramic copings or crowns shall be milled or polished with diamond instruments and with minimal pressure and water-cooling. The minimal thickness shall be 0,5 mm. Sharpe edges must be avoided.

Veneering: Copings shall be veneered with appropriate ceramics before cementing onto the abutment. The instructions for use of the ceramic manufacturers have to be considered.

Treatment of the titanium base and the ceramic abutment before cementing: Sandblasting of the contact surfaces with Al₂O₃, 50 µm, 2 bar and intensive cleaning of dust and grease. It is recommended to protect the connection part of the titanium base with an implant analog during handling.

Cementing: It is recommended to cement the ceramic abutment onto the titanium base with Panavia® F2.0 (Kuraray) or other equivalent cements. The instructions for use of the cements shall be followed carefully.

The titanium base shall be fixed onto an implant analog with the abutment screw. The head of the screw has to be covered with wax or resin. The mixed cement is applied onto the contact part of the titanium base. The abutment is pressed onto the ti-base. The final position is evaluated by slight rotation. The gap between abutment and base must be as small as possible. Remaining cement shall be removed immediately (*cementing instructions - only in German: www.cadstar.dental/Help Center/Technische Anleitungen/Klebeanleitung*).

Polishing: After hardening the remaining cement shall be removed with rotating silicon instruments. The cement inside the screw channel has to be removed carefully.

METAL DUST! Metal dust is harmful to your health. When milling and sandblasting, use a suction extraction system and a breathing mask.

Cleaning, disinfection and sterilization: The prosthetic components of our series are supplied in non-sterile condition. The instructions for cleaning and hygiene have to be considered, see detailed instruction manual GA_E_11-1-Serie in English (www.cadstar.dental/Help Center/Gebrauchsanweisungen).

Scan abutment: In the scanning procedure, the scan abutment is used to indicate the position of the implant. The size of the scan abutment shall be corresponding to the original implant system, implant diameter and titanium base series. The position of the beveled surface of the scan abutment corresponds to the position of the anti-rotation protection of the titanium base. The scan abutment is fixed on the implant analog with the abutment screw. After correct positioning and fixing, there is no gap visible between implant and scan abutment. Rotation of the scan abutment is impossible.

General information

Secondary effects: Allergies to the alloy, or contents of the alloy, or electrochemically based reactions may very rarely occur.

Reactions: In case of occlusal or approximal contact of different alloys, electrochemically based reactions may very rarely occur.

Warranty: 10 years on the mechanical stability of the titanium base, if it was processed according to the instruction for use.

If you have **any questions** about further **prosthetic components** of your series please contact our responsible member of staff:

Mr. Nikolaus Nehr / phone: +43 (0) 64 62 - 328 80 / e-mail: nikolaus.nehr@cadstar.dental