









Connections

✓ Supply possible
 ✗ Supply not possible
 ○ Supply only possible via model delivery

Connection	Platform	Ø	occlusally screwed bar/bridge	Individual titanium abutment	Zirconium abutment with titanium base
IB		Mini (2 mm)	✗	✓	✓
IB		Standard / Wide (3 mm)	✗	✓	✓
MU		Mini & Standard / Wide	✓	✗	✗

* only with original Argon K3®Pro titanium bases for Ø 2 mm stanchions ** only with original Argon K3®Pro titanium bases for Ø 3 mm stanchions

Astra Tech® Dental Implant System EV

IB		3,0 mm	✓	✓	✓
IB		3,6 mm	✓	✓	✓
IB		4,2 mm	✓	✓	✓
IB		4,8 mm	✓	✓	✓
IB		5,4 mm	✓	✓	✓

* only with original Astra Tech® prosthetics

Gingival height 1









Titanium Base EV Ø 3,0 mm (order no. 25924)
 Titanium Base EV Ø 3,6 mm (order no. 25926)
 Titanium Base EV Ø 4,2 mm (order no. 25928)
 Titanium Base EV Ø 4,8 mm (order no. 25930)
 Titanium Base EV Ø 5,4 mm (order no. 25932)

Gingival height 2

Titanium Base EV Ø 3,0 mm (order no. 25925)
 Titanium Base EV Ø 3,6 mm (order no. 25927)
 Titanium Base EV Ø 4,2 mm (order no. 25929)
 Titanium Base EV Ø 4,8 mm (order no. 25931)
 Titanium Base EV Ø 5,4 mm (order no. 25933)








Connections

✓ Supply possible
 ✗ Supply not possible
 ○ Supply only possible via model delivery





Connection	Platform	Ø	occlusally screwed bar/bridge	Individual titanium abutment	Zirconium abutment with titanium base
Astra Tech® Dental OsseoSpeed™					
IB		3,5 / 4,0 mm	✓	✓	✓
IB		4,5 / 5,0 mm	✓	✓	✓
MU		20°	○	✗	✗
MU		45°	○	✗	✗
Bego Semados® S-/RI-Line					
IB		3,25 / 3,75 mm	✓	✓	✓
IB		4,1 mm	✓	✓	✓
IB		4,5 mm	✓	✓	✓
IB		5,5 mm	○	✓	✓
<p>* only with original Bego Semados® CAD/CAM titanium base for S-/RI-Line (sub-system)</p> <p>S Ø 3.25/3.75 mm and RI Ø 3.75 mm (order no. 57285)</p> <p>S/RI Ø 4.1 mm (order no. 57287)</p> <p>S/RI Ø 4.5 mm (order no. 57289)</p> <p>S/RI Ø 5.5 mm (order no. 57291)</p>					

Connections

✓ Supply possible
 ✗ Supply not possible
 ○ Supply only possible via model delivery





Connection	Platform	Ø	occlusally screwed bar/bridge	Individual titanium abutment	Zirconium abutment with titanium base
Biomet 3i™ External					
IB		3,4 mm	✓	✓	✓
IB		4,1 mm	✓	✓	✓
IB		5,0 mm	✓	✓	✓
IB		6,0 mm	✓	✓	✓
MU		3,4 mm	✓	✗	✗
MU		4,1 / 5,0 mm	✓	✗	✗
MUD		MUD	○	✗	✗

Biomet 3i™ Certain®





IB		3,4 mm	✓	✓	✓
IB		4,1 mm	✓	✓	✓
IB		5,0 mm	✓	✓	✓
IB		6,0 mm	✓	✓	✓

Connections

✓ Supply possible
 ✗ Supply not possible
 ○ Supply only possible via model delivery

Connection	Platform	Ø	occlusally screwed bar/bridge	Individual titanium abutment	Zirconium abutment with titanium base
MU		3,4 mm	✓	✗	✗
MU		4,1 / 5,0 mm	✓	✗	✗
MU		6,0 mm	✓	✗	✗
MUD		MUD	○	✗	✗









Bredent Medical Sky®

IB		3,5 / 4,0 / 4,5 mm	✓	✓	✓
MU		3,5 / 4,0 / 4,5 mm (fast&fixed)	✓	✗	✗
IB		3,5 / 4,0 / 4,5 mm (uni.cone)	✓	✗	✗
MU		2,8 / 3,2 mm (miniSKY)	✓	✗	✗

* only with original Sky® uni.fit CAD abutment (order no. SKYUFCAD)

Connections

✓ Supply possible
 ✗ Supply not possible
 ○ Supply only possible via model delivery





Connection	Platform	Ø	occlusally screwed bar/bridge	Individual titanium abutment	Zirconium abutment with titanium base
Camlog Biotechnologies CAMLOG®					
IB		3,3 mm	✓	✓	✓
IB		3,8 mm	✓	✓	✓
IB		4,3 mm	✓	✓	✓
IB		5,0 / 6,0 mm	✓	✓	✓
MU		3,3 / 3,8 / 4,3 mm	✓	✗	✗
MU		5,0 / 6,0 mm	✓	✗	✗
MU		VSR 3,8 / 4,3 mm	○	✗	✗
MU		VSR 5,0 / 6,0 mm	○	✗	✗

Connections





✔ Supply possible
 ✘ Supply not possible
 ○ Supply only possible via model delivery

Connection	Platform	Ø	occlusally screwed bar/bridge	Individual titanium abutment	Zirconium abutment with titanium base
------------	----------	---	-------------------------------------	------------------------------------	---

Camlog Biotechnologies CONELOG®




IB		3,3 mm	✔	✔	✔
IB		3,8 mm	✔	✔	✔
IB		4,3 mm	✔	✔	✔
IB		5,0 mm	✔	✔	✔

Camlog Biotechnologies CONELOG®

MU		3,3 / 3,8 / 4,3 mm	✔	✘	✘
MU		5,0 mm	✔	✘	✘
MU		VSR 3,8 / 4,3 mm	○	✘	✘
MU		VSR 5,0 mm	○	✘	✘



Connections

✓ Supply possible
 ✗ Supply not possible
 ○ Supply only possible via model delivery

Connection	Platform	Ø	occlusally screwed bar/bridge	Individual titanium abutment	Zirconium abutment with titanium base
Dentaurum tioLogic®					
IB		S	✓	✓	✓
IB		M	✓	✓	✓
IB		L	✓	✓	✓




* only with original Wieland Zenotec® titanium base S tioLogic® (order no. 386-180-00)
 ** only with original Wieland Zenotec® titanium base M tioLogic® (order no. 386-280-00)
 *** only with original Wieland Zenotec® titanium base L tioLogic® (order no. 386-380-00)

Dentsply Friadent Ankylos® C/X

IB		3,5 mm	✓	✗	✓
MU		3,5 mm	✓	✗	✗

* only with original Dentsply Friadent Ankylos® TitaniumBase
 Ankylos® TitaniumBase C/ (not rotation locked / order no. 3102 2642)
 Ankylos® TitaniumBase /X (rotation locked / order no. 3102 2646)

Dentsply Friadent Frialit® & XiVE®

IB		3,4 mm	✓	✓	✓
IB		3,8 / 4,5 / 5,5 mm	✓	✓	✓
MU		3,4 / 3,8 / 4,5 / 5,5 mm	✓	✗	✗

Connections

✓ Supply possible
 ✗ Supply not possible
 ○ Supply only possible via model delivery

Connection	Platform	Ø	occlusally screwed bar/bridge	Individual titanium abutment	Zirconium abutment with titanium base
Henry Schein alphatech®					
IB		3,4 / 3,8 mm	✓	✗	✗
IB		4,3 / 5,0 mm	✓	✗	✗
IB		6,0 / 7,0 mm	✓	✗	✗
MU		3,4 / 3,8 / 4,3 / 5,0 / 6,0 / 7,0 mm	✓	✗	✗

Keystone PrimaConnex® & Genesis

IB		SD	✓	✓	✓
IB		RD	✓	✓	✓
IB		WD	✓	✓	✓










* only with original Keystone SD Titanium Base (order no. 700131)
 ** only with original Keystone RD Titanium Base (order no. 700132)
 *** only with original Keystone WD Titanium Base (order no. 700133)

Keystone Restore®

IB		SD	✓	✓	✗
IB		RD	✓	✓	✗
IB		WD	✓	✓	✗

Connections

✓ Supply possible
 ✗ Supply not possible
 ○ Supply only possible via model delivery

Connection	Platform	Ø	occlusally screwed bar/bridge	Individual titanium abutment	Zirconium abutment with titanium base
Medentis Medical templant®					
IB		3,75 / 4,1 / 4,8 mm	○	✗	✗
Medentis Medical ICX templant®					
IB		3,75 / 4,1 / 4,8 mm	✓	✓	✓
* only with original ICX templant® titanium base (order no. C-029-000002)					
Medical Instinct BoneTrust®					
MU		3,4 / 4,0 / 5,0 mm	○	✗	✗
Medentika MedentiBASE®					
MU		All Medentika MedentiBASE®	✓	✗	✗
Neoss Ltd. NEOSS™					
IB		3,5 / 4,0 / 4,5 / 5,0 / 5,5 mm	✓	✗	✓
MU		3,5 / 4,0 / 4,5 / 5,0 / 5,5 mm	✓	✗	✓
* only with original Neoss® NeoLink™ Titanium Mono titanium base (order no. 31133) ** only with original Neoss® Access NeoLink™ Titanium Multi titanium base (order no. 31226)					
Nobel Biocare Active® / (CC – Conical Connection) / Implant Direct InterActive™					
IB		NP	✓	✓	✓
IB		RP	✓	✓	✓
MU		NP / RP	✓	✗	✗

Connections





✓ Supply possible
 ✗ Supply not possible
 ○ Supply only possible via model delivery

Connection	Platform	Ø	occlusally screwed bar/bridge	Individual titanium abutment	Zirconium abutment with titanium base
------------	----------	---	-------------------------------------	------------------------------------	---

Nobel Biocare Brånemark®







IB		NP	✓	✓	✓
IB		RP	✓	✓	✓
IB		WP	✓	✓	✓
MU		NP / RP	✓	✗	✗
MU		WP	○	✗	✗

Nobel Biocare Replace Select® / Implant Direct Tri-Lobe

IB		NP	✓	✓	✓
IB		RP	✓	✓	✓
IB		WP	✓	✓	✓
MU		NP / RP	✓	✗	✗

Connections

✓ Supply possible
 ✗ Supply not possible
 ○ Supply only possible via model delivery

Connection	Platform	Ø	occlusally screwed bar/bridge	Individual titanium abutment	Zirconium abutment with titanium base
Paltop® Advanced®					
IB		NP (Narrow)	○	✗	✗
IB		SP (Standard)	○	✗	✗
MU		NP / SP	○	✗	✗
SIC invent AG SICase®					
IB		3,7 / 4,2 mm	✓	✗	✓
IB		4,7 / 5,0 mm	✓	✗	✓
MU		3,7 / 4,2 / 4,7 / 5,0 mm	✓	✗	✗





* only with original SIC Titanium titanium base for SICace® Ø 3.3 and Ø 4.2 mm (order no. 936190)
 ** only with original SIC Titanium titanium base for SICace® Ø 5.0 mm (order no. 936196)

Connections




✓ Supply possible
 ✗ Supply not possible
 ○ Supply only possible via model delivery

Connection	Platform	Ø	occlusally screwed bar/bridge	Individual titanium abutment	Zirconium abutment with titanium base
------------	----------	---	-------------------------------------	------------------------------------	---



Straumann BoneLevel®

IB		NC	✓	✓	✓
IB		RC	✓	✓	✓
MU		NC	✓	✗	✗
MU		RC	✓	✗	✗

Straumann synOcta® / Implant Direct Swish™

IB		NN	○	✗	✗
IB		RN	✓	✓	✓
IB		WN	✓	✓	✓





Prowital wi.tal®

IB		3,5 / 4,3 / 5,0 mm	✓	✗	✓
MU		3,5 / 4,3 / 5,0 mm	✓	✗	✗


* only with original Wi.tal® titanium base CAD/CAM Ø 3.5 (order no. 1.264.3500), Ø 4.3 (order no. 1.264.4300), Ø 5.0 (order no. 1.264.5000)

Connections


✓ Supply possible
 ✗ Supply not possible
 ○ Supply only possible via model delivery

Connection	Platform	Ø	occlusally screwed bar/bridge	Individual titanium abutment	Zirconium abutment with titanium base
Zimmer Dental Tapered Screw-Vent® / Implant Direct Legacy™ (Legacy 1-4)					
IB		3,5 mm	✓	✓	✓
IB		4,5 mm	✓	✓	✓
IB		5,7 mm	✓	✓	✓
MU		3,5 / 4,5 / 5,7 mm	✓	✗	✗

BTI Biotechnology Institute Interna®


IB		4,1 mm	○	✓	✗
----	---	--------	---	---	---

OT medical OT-F¹

IB		3,3 / 3,8 / 4,1 / 4,9 mm	✗	✗	✓
----	---	--------------------------	---	---	---

*only with original OT medical titanium base 3,30mm (order no. 01-5339002000CS), 3,80mm (order no. 01-5389002000CS), 4,10mm (order no. 01-5419002000CS), 4,90mm (order no. 01-5499002000CS)

OT medical OT-F² / OT-F³

IB		3,4 / 3,8 / 4,1 / 5,0 mm	✓	✓	✓
----	---	--------------------------	---	---	---

*only with original OT medical OT-F² titanium base 3,40mm (order no. 02-5349002010CS),

*only with original OT medical OT-F² / OT-F³ titanium base 3,80mm (order no. 02-5389002010CS), 4,10mm (order no. 02-5419002010CS), 5,00mm (order no. 02-5509002010CS)